

Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Typically used for solid-liquid separations in chemistry and biological laboratories, set-ups consist of a funnel, vacuum adapter and optional threaded connector for vials.

The vacuum adapters (CG-1049-J-series) are designed so that vacuum can be controlled with a PTFE stopcock to prevent filtrates from being pulled into the vacuum line. Stems on the filter funnels are 80mm long and prevent contamination of the vacuum adapter during the filtration process.

Part No.	Funnel Adapters with Vacuum Stopcocks	Price
CG-1049-J-14	Outer to 14/20 Inner	111.50 *
CG-1049-J-24	Outer to 24/25 Inner	114.90 *
CG-1049-J-29	Outer to 29/26 Inner	114.90 *
CG-1049-J-30	Vial, Outer to 20-400 GPI Thread	111.70 *

Part No.	Accessory Threaded Connectors	Price
CG-1318-23	20-400 to 20-400 GPI Thread	23.70 *
CG-1318-24	20-400 to 22-400 GPI Thread	24.30 *
CG-1318-25	20-400 to 24-400 GPI Thread	31.60 *

Disposable Filter Funnels with Polyethylene Frits

Available in Bulk Bags or Small-Quantity Trays


Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-10	18	10	Flared	100	290.00 *
OP-6602-12	40-60	10	Flared	100	318.75 *
OP-6602-13	40-60	40	Flared	100	321.25 *
OP-6602-14	110	10	Straight	50	243.75 *
OP-6602-16	300	20	Straight	30	273.75 *
OP-6602-18	500	20	Straight	20	256.25 *

Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-C10	18	10	Flared	24	84.60 *
OP-6602-C12	40	10	Flared	24	94.00 *
OP-6602-C14	110	10	Straight	12	74.00 *



Disposable Filter Funnels & Adapter System

Disposable Filter Funnels with Polyethylene and Diatomaceous Earth Frits

Available in Bulk Bags or Small-Quantity Trays

Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-10	18	500mg	Flared	50	245.00 *
OP-6603-12	40-60	1,500mg	Flared	50	261.25 *
OP-6603-13	150	5g	Straight	30	230.00 *

Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.


Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-C10	18	500mg	Flared	24	146.10 *
OP-6603-C12	40	500mg	Flared	24	156.50 *
OP-6603-C13	150	500mg	Straight	12	114.70 *



Cone-Shaped, Disposable Filter Funnels, Molded Polypropylene with Polyethylene Frits

Disposable Filter Funnels, Cone-Shaped, PP with PE Frits

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Case Qty.	Price
OP-6605-08	8	10	100	242.50 *
OP-6605-10	20	10	100	295.00 *
OP-6605-12	60	10	100	322.50 *



Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Purchase components separately.

Part No.	Funnel Adapters with Vacuum Stopcocks	Price
CG-1049-J-14	Outer to 14/20 Inner	111.50 *
CG-1049-J-24	Outer to 24/25 Inner	114.90 *
CG-1049-J-29	Outer to 29/26 Inner	114.90 *
CG-1049-J-30	Vial, Outer to 20-400 GPI Thread	111.70 *

Part No.	Accessory Threaded Connectors	Price
CG-1318-23	20-400 to 20-400 GPI Thread	23.70 *
CG-1318-24	20-400 to 22-400 GPI Thread	24.30 *
CG-1318-25	20-400 to 24-400 GPI Thread	31.60 *

