



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number 220311

PCN Title CG-2077-H Bearings

Notification Date 3/11/2022

Implementation Date 3/11/2022

Part Description Change Compression Cap to Black Peek Material

Affected Part Number(s) CG-2077-H-01 (Item CG-2077-H-20), CG-2077-H-02 (Item CG-2077-H-20)
CG-2077-H-05 (Item CG-2077-H-21), CG-2077-H-06 (Item CG-2077-H-22)
CG-2077-H-07 (Item CG-2077-H-23)

Change Description Changed Compression Cap from Black Nylon / Acetal to Black Peek

Reason for Change Changed for a better Chemical Compatibility Material.

Performance Impact Material Change will increase the Resistance of Deterioration in Corrosive Environments.

PCN Correspondence QC@cglifesciences.com

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.