



3800 N. Mill Road  
Vineland NJ 08360  
Phone: 800.843.1794  
Fax: 800.922.4361  
[www.cglifesciences.com](http://www.cglifesciences.com)

## PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

**PCN Number** 240215

**PCN Title** CG-309-B, CG-309-U, CG-309-Z

**Notification Date** 2/15/2024

**Implementation Date** 2/16/2024

**Part Description** FFKM Perfluoro O'Rings Black, White, Class VI 486 FFKM

**Affected Part Number(s)** CG-309-Series, CG-309-B Series, CG-309-U Series, CG-309-Z Series

**Change Description** CG-309-FFKM Perfluoroelastomer has been discontinued by the manufacturer. This FFKM material is being replaced by CG-309-U White FFKM Perfluoro compound and CG-309-B Black FFKM Perfluoro compounds.

**Reason for Change** CG-309-Series Perfluoroelastomer has been discontinued by the manufacturer.

**Performance Impact** No Impact.

**PCN Correspondence** [QC@cglifesciences.com](mailto:QC@cglifesciences.com)

# PRODUCT CHANGE NOTIFICATION

## FORM - FIT - FUNCTION DEFINITIONS

**FORM:** Refers to such characteristics as external dimensions, weight and size of a part or assembly.

**FIT:** Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

**FUNCTION:** Is a criteria that is met when the part performs its stated purpose effectively and reliably.