

Mix 'n' Match
**Buy Two
 Get Third
 FREE!**
 Promo Code **ACCUPET23**

AccuPet Pro™ Micropipettes

- digital accuracy
- continuously adjustable
- research grade
- designed for comfort

rubberized grips

universal tip fitting shafts

autoclavable lower portion

DirectVue volume indicators

Comfort

Function

Affordability

AccuPet - The Standard for Microliter Measurement

Part No.	Volume Range	Accuracy %	Precision %	Price
CLS-3489E-01	0.1-2.5µL	5	1.2	198.00
CLS-3489E-02	0.5-10µL	1	0.8	198.00
CLS-3489E-03	2-20µL	1	0.4	198.00
CLS-3489E-04	5-50µL	0.6	0.5	198.00
CLS-3489E-05	10-100µL	0.8	0.2	198.00
CLS-3489E-06	20-200µL	0.8	0.2	198.00
CLS-3489E-07	100-1,000µL	0.8	0.2	198.00
CLS-3489E-08	500-5,000µL	0.6	0.2	198.00
CLS-3489E-09	500-10,000µL	0.6	0.2	198.00

Note: Accuracy and precision value given at maximum volume of pipettor. No tip ejector on the 10mL model.



Place your order today!

www.cglifesciences.com
 customer-service@cglifesciences.com
 technical-service@cglifesciences.com

T 1.800.843.1794

F 1.800.922.4361

Printed in USA. CG-23074IHpdf0923

* Offer valid from 1/4/2023 through 12/31/2023. Must purchase at list price and specify promo code **ACCUPET23** when placing the order. Valid on US and Canadian purchases only. May not be combined with any other offer or discount. This offer

is non-transferable and void where prohibited by law or company policy. Offer is subject to Chemglass Life Sciences' standard terms, conditions and contracts Offer not valid through distributors.