

# Advanced Cloud Control for Reactors

## NEW Upgraded Features

Get realtime monitoring and notifications!



### Chemglass / Velp Wi-Fi Enabled, Digital Motor Controller (CG-2033-V-75)

- Safety "Lock-Out" function
- Real time torque and speed graph
- Built-in timer and PT 100 temp monitor
- Nemko mark indicates CSA and UL compliance
- Cloud enabled - Ermes realtime monitoring and notifications
- Upgrade - compatible with existing CG-2033-B-25 motor
- USB interface for PC software data log, "Off / On" monitoring

Get all the details here:

[www.cglifesciences.com](http://www.cglifesciences.com)



### Sample Screens



**Ermes Cloud One-Year Subscription**  
Part No. CG-1997-W-601



tech@cglifesciences.com  
T 1.800.843.1794 F 1.800.922.4361  
[www.cglifesciences.com](http://www.cglifesciences.com)

# Advanced Cloud Control for Reactors

OHS USB Software




VS.

Velp Ermes Cloud

Compare features of standard OHS USB software with the NEW Cloud-enabled Controller below.


## OHS USB Software (CG-2033-V-600)

- Plug and play software for data analysis
- Local data logging, easy to read real-time graph
- Record up to four instruments at once
- Export data to .csv or load onto OHS



## NEW Velp Ermes Cloud (CG-1997-W-601)

- Easy-to-use Wi-Fi communication
- Always connected to historical data and live monitoring via cloud
- Immediate notifications and alerts on mobile devices or computer
- OTA software updates for added security and features
- Unlimited instrument monitoring
- Create reports instantly in many formats such as .pdf, .txt, .xls, LIMS, etc.
- Separate instruments by lab, hood, benchtop, etc.
- Safely manage data with maximum level of encryption, following the highest cyber-security standards
- User management with permission settings
- Monitor and manage multiple instruments and consumables 24 / 7 from your pc, tablet or smartphone



## Compact Design



WiFi Enabled Controller  
Part No.  
CG-2033-V-75



## Sample Screens



tech@cglifesciences.com  
T 1.800.843.1794 F 1.800.922.4361  
www.cglifesciences.com

CG-19082IH-pdf1019