

# PRODUCT CHANGE NOTIFICATION

**PCN Number** 010220

**PCN Title** Solvent Still Head

**Notification Date** 1/22/2020

**Implementation Date** 1/22/2020

**Part Description** 125ml and 500ml Solvent Still Head

**Affected Part Number(s)** AF-0700-01, 125ml Solvent Still Head and AF-0700-02, 500ml Solvent Still Head

**Change Description** Changed angle of 22 m/m x 1.8 Wall from 30° to 45°.

**Reason for Change** Customers complain they cannot wrap the tube because there was not enough space between the tube and the body.

**Performance Impact** NONE

**Distinguishable Characteristics** Tube is now on a 45° Angle from Body.

**PCN Correspondence** [QC@cglifesciences.com](mailto:QC@cglifesciences.com)

# PRODUCT CHANGE NOTIFICATION

## FORM - FIT - FUNCTION DEFINITIONS

**FORM:** Refers to such characteristics as external dimensions, weight and size of a part or assembly.

**FIT:** Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

**FUNCTION:** Is a criteria that is met when the part performs its stated purpose effectively and reliably.