

Pressure Relief, Storage and Sample

VIALS



Pressure Relief Reaction Vials are perfect for use in heated reactions where internal pressure created by the reaction could exceed 150 psig.

Pressure Relief Reaction Vials in Cell Packs

Supplied with special open top caps to prevent internal pressure from exceeding 150 psig. Made from 33 expansion borosilicate glass for high resistance to thermal shock and maximum chemical resistance. Conform to ASTM E-438 Type 1, Class A, USP Type 1 and ISO 3585. Packaged in SBS, coated, white trays complete with SBS partitions. Boxes are marked with alpha-numeric coordinates for fast and easy vial identification. Packed 100 per box.

Part Number	Vial Capacity	Vial Dimensions	Thread Size	Pkg. Price
CG-4912-01	1 dram	15 x 45	13-425	93.05 ▲
CG-4912-02	2 dram	17 x 60	15-425	101.05 ▲
CG-4912-05	20mL	28 x 57	24-414	109.15 ▲
CG-4912-06	40mL	28 x 95	24-414	144.80 ▲



CG-4912 Pressure Relief

Solid Top Storage Vials in Cell Packs

For storage applications, these vials are packaged in SBS, coated, white trays complete with SBS partitions. Boxes are marked with alpha-numeric coordinates for fast and easy vial identification. Caps are supplied attached and all sizes have the TFE-lined, silicone, SURE-LINK septa. The 4mL is a 1 dram, the 8mL is a 2 dram. Packed 100 per box.

Part Number	Vial Capacity	Vial Dimensions	Thread Size	Cap Color	Pkg. Price
CG-4908-04	4mL (1 dram)	15 x 45	13-425	Green	63.50 ▲
CG-4908-05	8mL (2 dram)	17 x 60	15-425	Green	69.20 ▲
CG-4908-03	20mL	28 x 57	24-400	White	78.90 ▲
CG-4909-07	40mL	28 x 95	24-400	White	115.65 ▲



CG-4908 / CG-4909 Solid Tops

Open Top Storage Vials with Septa in Cell Packs

For reaction or storage applications, these vials are packaged in SBS, coated, white trays complete with SBS partitions. Boxes are marked with alpha-numeric coordinates for fast and easy vial identification. Caps are supplied attached and all sizes have the TFE-lined, silicone, SURE-LINK septa. Packed 100 per box.

Part Number	Vial Capacity	Vial Dimensions	Thread Size	Cap Color	Pkg. Price
CG-4909-04	4mL	15 x 45	13-425	Green	56.10 ▲
CG-4909-03	8mL	17 x 60	15-425	Green	63.50 ▲
CG-4908-01	20mL	28 x 57	24-400	Green	71.55 ▲
CG-4909-05	40mL	28 x 95	24-400	Green	108.35 ▲



CG-4908 / CG-4909 Open Tops

▲ Subject to Vials / Containers Discount

Clear Sample Vials with Solid Top Caps, LAB-PAC

Sample vials are packaged in convenient LAB-PACS of 144 each. Vials are made from Type 1, borosilicate glass. Full height partitions minimize breakage. Solid, PTFE-lined caps are attached.

Part Number	Vial Capacity	Vial Dimensions	Thread Size	Pkg. Price
CG-4900-01	4mL (1 dram)	15 x 45	13-425	136.60
CG-4900-02	8mL (2 dram)	17 x 60	15-425	144.60
CG-4900-03	16mL (4 dram)	21 x 70	18-400	156.75
CG-4900-04	20mL	28 x 57	24-400	230.25
CG-4900-05	30mL	28 x 70	24-400	236.85
CG-4900-06	40mL	28 x 95	24-400	257.00

Replacement Open Top Screw Caps with PTFE-Faced Silicone Septa

Part Number	Thread Size	Color	Pkg. Qty.	Pkg. Price
CG-4910-01	13-425	black	144	48.60
CG-4910-02	15-425	black	144	68.45
CG-4910-03	18-400	black	144	90.60
CG-4910-04	24-400	green	144	145.95



CG-4900 Solid Tops

Chemglass Stirring Magnetic Hotplates and Accessories

High Performance • Feature-Packed • Affordably Priced

Pie-Block Safety System

CG-1991-P-01
\$408.55



Analog with 3-Year Warranty

CG-1996-V10
\$545.80



Digital with Probe and 5-Year Warranty

CG-1997-V-10 \$745.00



T 1.800.843.1794

F 1.800.922.4361

customer-service@cglifesciences.com

www.cglifesciences.com

▲ Subject to Vials / Containers Discount CG-17079BC7.5M0817