

Individual Rolls of Colorful Labeling Tape for Glass or Plastic

Choose from
Single Rolls
OR
**Rainbow
Packs**

14-Roll RAINBOW Packs



RAINBOW Packs include
two rolls each of seven
colors: white, yellow, green,
red, orange, blue and
lavender.



- 3/4- and 1-inch rolls
- adheres to any dry surface
- removes cleanly and easily
- durable, water-proof, oil and acid-resistant
- withstands temperatures from -23 to 121°C
- can be written on with permanent marker, pencil or pen
- 500" rolls with 1" cores

Single Rolls in **31 colors** listed on back ►



T 1.800.843.1794

F 1.800.922.4361

www.cglifesciences.com

customer-service@cglifesciences.com



**Single
Rolls or
Rainbow
Packs**



Rainbow Packs
are listed at the
bottom of the
page.

Choose Single Roll Tape Colors		No. of Rolls	3/4" x 500" Roll		No. of Rolls	1" x 500" Roll	
			Item Number	Each Price		Item Number	Each Price
	white		CG-4035-21	\$6.55*		CG-4035-61	\$8.05*
	yellow		CG-4035-22	\$6.55*		CG-4035-62	\$8.05*
	green		CG-4035-23	\$6.55*		CG-4035-63	\$8.05*
	red		CG-4035-24	\$6.55*		CG-4035-64	\$8.05*
	orange		CG-4035-25	\$6.55*		CG-4035-65	\$8.05*
	blue		CG-4035-26	\$6.55*		CG-4035-66	\$8.05*
	pink		CG-4035-27	\$6.55*		CG-4035-67	\$8.05*
	tan		CG-4035-28	\$6.55*		CG-4035-68	\$8.05*
	salmon		CG-4035-29	\$6.55*		CG-4035-69	\$8.05*
	copper		CG-4035-30	\$6.55*		CG-4035-70	\$8.05*
	rose		CG-4035-31	\$6.55*		CG-4035-71	\$8.05*
	violet		CG-4035-32	\$6.55*		CG-4035-72	\$8.05*
	lavender		CG-4035-33	\$6.55*		CG-4035-73	\$8.05*
	gray		CG-4035-34	\$6.55*		CG-4035-74	\$8.05*
	lime		CG-4035-35	\$6.55*		CG-4035-75	\$8.05*
	aqua		CG-4035-36	\$6.55*		CG-4035-76	\$8.05*
	chartreuse		CG-4035-37	\$6.55*		CG-4035-77	\$8.05*
	gold		CG-4035-38	\$6.55*		CG-4035-78	\$8.05*
	silver		CG-4035-39	\$6.55*		CG-4035-79	\$8.05*
	magenta		CG-4035-40	\$6.55*		CG-4035-80	\$8.05*
	dark red		CG-4035-41	\$6.55*		CG-4035-81	\$8.05*
	maroon		CG-4035-42	\$6.55*		CG-4035-82	\$8.05*
	dark brown		CG-4035-43	\$6.55*		CG-4035-83	\$8.05*
	dark blue		CG-4035-44	\$6.55*		CG-4035-84	\$8.05*
	dark green		CG-4035-45	\$6.55*		CG-4035-85	\$8.05*
	black		CG-4035-46	\$6.55*		CG-4035-86	\$8.05*
	fluorescent red		CG-4035-47	\$6.55*		CG-4035-87	\$8.05*
	fluorescent orange		CG-4035-48	\$6.55*		CG-4035-88	\$8.05*
	fluorescent green		CG-4035-49	\$6.55*		CG-4035-89	\$8.05*
	fluorescent chartreuse		CG-4035-50	\$6.55*		CG-4035-90	\$8.05*
	fluorescent pink		CG-4035-51	\$6.55*		CG-4035-91	\$8.05*

14 Roll Rainbow Packs (Two Rolls Each of White, Yellow, Green, Red, Orange, Blue and Lavender)

Choose 3/4" or 1" Rolls ►	CG-4035-R075	87.65*	CG-4035-R100	108.75*
---------------------------	--------------	--------	--------------	---------

* Not subject to discount.

©2021 Chemglass Life Sciences.
Printed in USA. Prices subject to
change without notice.

CG-21017IHpdf0221



T 1.800.843.1794

F 1.800.922.4361

www.cglifesciences.com

customer-service@cglifesciences.com