



3800 N. Mill Road  
Vineland NJ 08360  
Phone: 800.843.1794  
Fax: 800.922.4361  
www.cglifesciences.com

## PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

**PCN Number** 200214

**PCN Title** CG-1949-F104

**Notification Date** 2/14/2020

**Implementation Date** 2/14/2020

**Part Description** 100mm Filter Base, Side Outlet

**Affected Part Number(s)** CG-1949-F104

**Change Description** Would like to lift the o'ring seat by 1/32 (.031")

**Reason for Change** It was discovered in the original design the top flange bottoms out on the pinch ring before it makes contact with the o'ring causing a leak. By lifting the o'ring seat, it will allow the o'ring to seal properly.

**Performance Impact** Improved seal.

**PCN Correspondence** [QC@cglifesciences.com](mailto:QC@cglifesciences.com)