

Pressure Relief Reaction Vials



NEW!
Smaller packs
for easier
handling and
storage!

Packs of 24 and
48 vials each

Trays marked with
alpha-numerics



*See all sizes and pack
configurations on the other side.*

Convenient Shelf Packs of 24 & 48 Vials *For Easier Transport, Storage and Use in the Lab*

Chemglass Pressure Relief Reaction Vials are designed for optimal performance and safety in heated reactions.

- special open top caps prevent vial pressure from exceeding 150 psig
- borosilicate glass reaction vials for heat- and chemical-resistance*
- packaged in SBS-coated, white trays with SBS partitions
- trays are marked with alpha-numeric coordinates for fast and easy vial identification

* Made from 33 expansion borosilicate glass conforming to ASTM E-438 Type 1-Class A, USP Type 1 and ISO 3585

Now available in NEW, smaller quantity packs for easier transport, storage and use.



1.800.843.1794

cs@cglifesciences.com

www.cglifesciences.com

Pressure Relief Reaction Vials



Packaged to Fit Your Needs

Designed for Optimal Performance and Safety in Heated Reactions

NEW Convenient Shelf Packs of 24 & 48 Vials

Smaller Quantities For Easier Transport, Storage and Use in the Lab

2 Dram Reaction Vials in Convenient Shelf Packs

Part No.	Vial Capacity	Vial Dims. (mm)	Thread Size	Qty.	2021 Price
CG-4912-242	2 dram	17 x 60	15-425	24	\$31.50 ▲
CG-4912-482	2 dram	17 x 60	15-425	48	\$62.50 ▲

20mL Reaction Vials in Convenient Shelf Packs

Part No.	Vial Capacity	Vial Dims. (mm)	Thread Size	Qty.	2021 Price
CG-4912-245	20mL	28 x 57	24-414	24	\$33.75 ▲
CG-4912-485	20mL	28 x 57	24-414	48	\$67.50 ▲

Larger Quantity 100-Packs

For Labs Needing Larger Quantities of Vials per Tray

Reaction Vials in Trays of 100 Vials Each

Part No.	Vial Capacity	Vial Dims. (mm)	Thread Size	Qty.	2021 Price
CG-4912-01	1 dram	15 x 45	13-425	100	\$100.70 ▲
CG-4912-02	2 dram	17 x 60	15-425	100	\$109.35 ▲
CG-4912-05	20mL	28 x 57	24-414	100	\$118.15 ▲
CG-4912-06	40mL	28 x 95	24-414	100	\$156.65 ▲

Purchase caps and vials complete. Not available separately.

▲ Subject to vials / containers discount.

Prices subject to change without notice. Printed in USA. CG-20092BC5M1220



Convenient Shelf Packs of 2 dram vials



Convenient Shelf Packs of 20mL vials

Larger Quantity 100-pack



1.800.843.1794

cs@cglifesciences.com

www.cglifesciences.com