



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number 210223

PCN Title CG-1413-SERIES

Notification Date 2/23/2021

Implementation Date 2/23/2021

Part Description Buchner Funnel, Table Tops

Affected Part Number(s)

- CG-1413-05**, Buchner Funnel, Table Top, 457mm, Medium Porosity Fixed Plate
- CG-1413-06**, Buchner Funnel, Table Top, 457mm, Coarse Porosity Fixed Plate
- CG-1413-07**, Buchner Funnel, Table Top, 457mm, Medium Porosity Removable Plate
- CG-1413-08**, Buchner Funnel, Table Top, 457mm, Coarse Porosity Removable Plate
- CG-1413-19**, Buchner Funnel, Table Top, 457mm, Perforated Fixed Plate
- CG-1413-20**, Buchner Funnel, Table Top, 457mm, Perforated Removable Plate

Change Description Product changed in height and shape. Elimination of all external structural welds.

OLD HEIGHT	NEW HEIGHT
11.50 in.	11.25 in.

Reason for Change This is an aesthetic and structural change. Product is 46% lighter in weight than previous model. Higher structural integrity.

Performance Impact Product is lighter in weight with higher structural integrity. Stackable design for easy storage.

PCN Correspondance QC@cglifesciences.com

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.