



3800 N. Mill Road  
Vineland NJ 08360  
Phone: 800.843.1794  
Fax: 800.922.4361  
[www.cglifesciences.com](http://www.cglifesciences.com)

## PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

<b>PCN Number</b>	<u>200730</u>
<b>PCN Title</b>	<u>CLS-1481-22</u>
<b>Notification Date</b>	<u>7/30/2020</u>
<b>Implementation Date</b>	<u>8/3/2020</u>
<b>Part Description</b>	<u>Cap, GL-45, Green, Open Top, ( 1-5/16" Dia. Hole) Polypropylene</u>
<b>Affected Part Number(s)</b>	<u>N/A</u>
<b>Change Description</b>	<u>Chemglass to manufacture cap. Color to be slightly different and Chemglass Bio Sciences Logo on cap.</u>
<b>Reason for Change</b>	<u>Improved lead time and part tolerance due to one step molding process and eliminates the machining process the previous part needed.</u>
<b>Performance Impact</b>	<u>N/A</u>
<b>PCN Correspondence</b>	<u><a href="mailto:QC@cglifesciences.com">QC@cglifesciences.com</a></u>

# PRODUCT CHANGE NOTIFICATION

## FORM - FIT - FUNCTION DEFINITIONS

**FORM:** Refers to such characteristics as external dimensions, weight and size of a part or assembly.

**FIT:** Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

**FUNCTION:** Is a criteria that is met when the part performs its stated purpose effectively and reliably.