



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number 210520

PCN Title Hour Meters to be added to CLS Equipment

Notification Date 5/20/2021

Implementation Date 5/20/2021

Part Description Digital Hour Meter

Affected Part Number(s) CLS-4100-04, Stirrer, Magnetic, Five Position, Compact
CLS-4100-A4, Stirrer, Magnetic, Five Position, Compact, Alarm
CLS-4100-06, Stirrer, Magnetic, Five Position
CLS-4100-A6, Stirrer, Magnetic, Five Position, Alarm
CLS-4100-09, Stirrer, Magnetic, Nine Position
CLS-4100-A9, Stirrer, Magnetic, Nine Position, Alarm
CLS-4021-100, Orbital Shaker, 120v, Versa-Orb, Base Only, Digital

Change Description Addition of hour meters to current CLS items.
NOTE: THESE AFFECTED PRODUCTS WILL UNDERGO A PHASE IN
PROCESS WITH THIS CHANGE.

Reason for Change To record and ultimately develop preventative maintenance
recommendations.

Performance Impact None

PCN Correspondence QC@cglifesciences.com

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.